


# Envoi Pen Injector

## DEVICE SPECIFICATIONS



152 mm total length (including cap)  
 Max elongation during dosing: +33 mm (delta)



Platform assumes cartridge size 3 ml; customizable for 1.5 ml



Acceptable stopper lengths: between 8 mm and 11 mm

17 mm diameter  
 3 ml inner diameter 9.70 mm



Standard platform color scheme



Custom platform color scheme on request

Push button color customizable to signal the different drug types

Intuitive, easy-scroll, dose-setting unit for variable or fixed dosing

Dosing print customizable

Light grey pen and cap

Optional translucent or full color cartridge housing with inspection window

### Bringing Possibilities To Life

Disclaimer: This product has not obtained a CE marking applicable to the health, safety and environmental regulations of the European Union, and has not been reviewed or cleared for use by the U.S. Food and Drug Administration or comparable health regulatory authorities in other jurisdictions.

© 2025 Phillips Medisize - a Molex Company. All Rights Reserved. PM-PH-1482041304-AA