


Envoi Pen Injector

GLP-1 – DEVICE SPECIFICATIONS



152 mm total length
(including cap)
Max elongation during
dosing: +33 mm (delta)



Platform assumes
cartridge size 3 ml;
customizable for 1.5 ml



Acceptable stopper lengths:
between 8 mm and 11 mm

17 mm diameter

3 ml inner diameter 9.70 mm



Standard platform color
scheme for GLP-1

The push button in a range of
colors to match the GLP-1 products
for diabetes care

Intuitive, easy-scroll, dose-setting
unit for seamless adjustment of
variable or fixed dosing

Dosing print matching
the GLP-1 requirements

Light grey pen housing and cap.
Customizable color available
upon request

The cartridge holder platform
assumes it comes translucent with
inspection windows
Customizable color available
upon request.

Bringing Possibilities To Life

Disclaimer: This product has not obtained a CE marking applicable to the health, safety and environmental regulations of the European Union, and has not been reviewed or cleared for use by the U.S. Food and Drug Administration or comparable health regulatory authorities in other jurisdictions.

© 2025 Phillips Medisize - a Molex Company. All Rights Reserved. PM-PH-1482041304-AA